

+ Improving incisional outcomes, for better lives following vascular surgery

PICO[◇] sNPWT has shown to significantly reduce the odds of developing a surgical site infection (SSI) by 78%¹ and reducing the incidence of seroma,² while saving cost^{3,4} following peripheral vascular surgery.*

Smith+Nephew



PICO[◇]
Single Use Negative Pressure
Wound Therapy System

www.smith-nephew.com/pico

Helping you get **CLOSER TO ZERO[◇]**
surgical site complications¹

*Compared to care with standard dressings; p=0.03; meta-analysis of 2 studies (OR: 0.22)

DISCOVER BETTER OUTCOMES HERE




How anatomy can dictate incisional outcomes

Open incisions in the inguinal region for vascular surgery can:

- **disrupt the lymphatic structures**, which can continue to leak for a long time after surgery^{5,6} be subject to increased lateral tension that cause stress on the suture lines following closure⁵
- **risk colonisation stemming from the skin flora**, proximity to external genitalia and relative moisture of the groin⁷

making them at **elevated risk of developing a post-operative SSI**.^{5,6}



Are your patients at elevated risk?



The risk of developing a post-operative SSC depends on the type of surgery and patient risk factors^{12,13}

The presence of just **1 major risk factor** or **2** or more moderate risk factors, places patients at elevated risk of an SSC and means you should consider **PICO^o sNPWT**.¹²



Category

Patient-related risk factor

Procedural-related risk factor



Major risk factor
presence of 1 = high risk of surgical site complication

! BMI $\geq 40\text{kg/m}^2$ or $\leq 18\text{kg/m}^2$

! Extended surgery*

! Uncontrolled insulin dependent diabetes mellitus

! Emergency surgery

! Renal dialysis

! Hypothermia



Moderate risk factor
presence of 2 \geq high risk of surgical site complication

! ASA physical status $>II$

! Anaemia / blood transfusion

! Age < 1 year or > 75 years

! High wound tension after closure

! BMI 30–39.9 kg/m^2

! Dual antiplatelet treatment

! Immunosuppression

! Suboptimal timing or omission of prophylactic antibiotics

! Smoking (current)

! Tissue trauma / large area of dissection / large area of undermining

! Pre-existing infection at a body site remote from operative site

Table adapted from World Union of Wound Healing societies Consensus, 2016. The risk factors represented in this table are examples only and not an exhaustive list¹² *Defined as $>T$ (hours) which is dependent on the type of surgical procedure, and is the 75th centile of duration of surgery for a particular procedure, e.g. coronary artery bypass graft has a T of 5 hours and caesarean section has a T of 1 hour.



Reducing the risk of complications with vascular surgery

Ask for

Evidence in focus
publication summary



In a randomised controlled trial (RCT) of 178 patients, the prophylactic use of **PICO^o sNPWT** **significantly reduced the incidence of SSIs by 64%** following vascular surgery in the inguinal region.*¹⁴

FOLLOWING VASCULAR SURGERY IN THE INGUINAL REGION



**64% relative reduction
in SSC incidence**

with PICO sNPWT versus standard dressings for pooled results (ASEPSIS criteria)

* Compared to care with standard dressings; p=0.02; NOTE: N=178 (pooled unilateral and bilateral incisions)



HOW ANATOMY DICTATES OUTCOMES



BACK AT HOME, NOT IN HOSPITAL



THE PICO SYSTEM



ORDERING



REFERENCES

PICO^o sNPWT FOR VASCULAR SURGERY | 04

Back at home, not back in hospital, after vascular surgery

Ask for

Evidence in focus
publication summary

PICO^o sNPWT helped to **significantly reduce the incidence of SSCs**, including an **82% relative reduction in the incidence of seroma**,* in patients undergoing vascular surgery in the inguinal region.²

Mean hospital length of stay[†] and time to resolution of wound complications[‡] **were shorter with PICO sNPWT** than with standard dressings for readmitted patients, **which contributed to cost savings**.^{||2}

FOLLOWING VASCULAR SURGERY IN THE INGUINAL REGION

**57% relative reduction
in wound complications**

with PICO sNPWT versus care with standard dressings

**82% relative reduction
in the incidence of seroma**

with PICO sNPWT versus care with standard dressings

* Compared to care with standard dressings † PICO 52 days; standard dressings 96 days

‡ PICO sNPWT (3 patients, 2.83 days) v standard dressings (6 patients, 5.67 days, p=0.465) || compared with standard dressings




Lowering the incidence of SSIs after vascular surgery

Ask for

Evidence in focus
publication summary

When all costs were considered **PICO[®] sNPWT** was determined to be cost effective* as it was demonstrated to help **significantly reduce the incidence of SSIs by 60%** when used prophylactically following vascular surgery in the inguinal region.^{†4}

FOLLOWING VASCULAR SURGERY IN THE INGUINAL REGION



**60% relative reduction
in the incidence of SSIs**

with PICO sNPWT versus standard dressings at 90 days postoperatively

* estimated incremental cost-effectiveness ratio, €1,853 per SSI avoided † compared with standard dressings, p=0.015



HOW ANATOMY DICTATES OUTCOMES



BACK AT HOME, NOT IN HOSPITAL



THE PICO SYSTEM



ORDERING



REFERENCES

PICO[®] sNPWT FOR VASCULAR SURGERY | 06

Reducing risks and costs after femoral endarterectomy

Ask for

Evidence in focus
publication summary

Following femoral endarterectomy, the prophylactic use of PICO[®] sNPWT helped **significantly reduce the rate of SSCs by 63%*** including dehiscence[†] and seroma.^{‡15}

PICO therapy was also calculated to be cost saving, **reducing the average per patient treatment cost by \$757.00 CDN***.^{§15}

“Clinicians found the dressing and device easy to apply and operate[¶]”

Taken directly from Wikkeling et al., 2021¹⁵

FOLLOWING FEMORAL ENDARTERECTOMY

**63% relative reduction
in the incidence of SSCs**

with PICO sNPWT versus care with standard dressings*

**\$757.00* CDN
estimated
cost saving per patient**

with PICO sNPWT versus care with standard dressings

* Conversion as of July 2022

*Compared with standard care; N=108; 50% v 18.18% p=0.0011 † Compared with standard care; 32.8% v 9.1% ‡ Compared with standard care; 10.9% v 4.5%

§ Compared with standard care ¶ All patients were operated on by the same two surgeons



HOW ANATOMY DICTATES OUTCOMES



BACK AT HOME, NOT IN HOSPITAL



THE PICO SYSTEM



ORDERING



REFERENCES

PICO[®] sNPWT FOR VASCULAR SURGERY | 07

The PICO[◇] 7 System

Completely portable and clinically effective in the treatment of surgical,¹ chronic^{16,17} and acute^{18,19} wounds.

PICO[◇] 7

Single Use Negative Pressure Wound Therapy System



Features:

Improved device performance*

- The PICO 7 pump has a significantly higher maximum leak rate tolerance than the original PICO pump^{†20}

Improved ease-of-use

- New user interface with a 'dressing full' indicator, optimising dressing changes²¹
- Area to write start date of therapy, helping with healthcare protocols²²

Designed to improve patient quality of life

- Now even quieter pump than before²³
- New transparent belt clip designed to enable greater portability

Increased flexibility

- New multipacks of five dressings now available, allowing therapy to be tailored to patients' clinical needs

* Compared with PICO † p < 0.001



HOW ANATOMY DICTATES OUTCOMES



BACK AT HOME, NOT IN HOSPITAL



THE PICO SYSTEM



ORDERING









REFERENCES

PICO[◇] sNPWT FOR VASCULAR SURGERY | 08

Product ordering codes

The PICO[◇] sNPWT portfolio is compatible^{16,24} with ACTICOAT[◇] FLEX Antimicrobial Barrier Dressing, our silver-coated antimicrobial wound contact layer.*

		PICO 7 System		PICO 14 System	Multipack with	PICO 7Y device
		+ 1 dressing	+ 2 dressings	+ 2 dressings	5 dressings	+ 2 dressings
Dressing sizes		Code	Code	Code	Code	Code
	Multisite small 15cm x 20cm	66802010	66802000	66802040	66802020	–
	Multisite large 20cm x 25cm	66802011	66802001	66802041	66802021	66802031
	10cm x 20cm	66802012	66802002	66802042	66802022	–
	10cm x 30cm	66802013	66802003		66802023	–
	10cm x 40cm	66802014	66802004		66802024	–
	15cm x 15cm	66802015	66802005	66802045	66802025	–
	15cm x 20cm	66802016	66802006	66802046	66802026	–
	15cm x 30cm	66802017	66802007	66802047	66802027	–
	20cm x 20cm	66802018	66802008	66802048	66802028	–
	25cm x 25cm	66802019	66802009	66802049	66802029	–
Consumables		Code				
	Foam dressing filler		10cm x 12.5cm	66801021		

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

*ACTICOAT FLEX 3 and FLEX 7 are only approved for use with NPWT for up to 3 days

Advanced Wound Management
Smith & Nephew Inc.
2280 Argentia Rd.
Mississauga, ON L5N6H8
T 1-800-463-7439
F 1-800-671-9140

www.smith-nephew.com

◇ Trademark of Smith+Nephew
All Trademarks acknowledged
© January 2022 Smith+Nephew
AWM-AWD-33653 | CA48480 08/22

Helping you get **CLOSER TO ZERO** surgical site complications¹

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

24/7 Negative Pressure Wound Therapy Hotline
1-800-463-7439

References: 1. Saunders C, Nherera LM, Horner A, Trueman P. Single-use negative-pressure wound therapy versus conventional dressings for closed surgical incisions: systematic literature review and meta-analysis. *BJs open*. 2021 Jan;5(1):zraa003. 2. Fleming CA, Kuteva M, O'Hanlon K, O'Brien G, McGreal G. Routine use of PICO dressings may reduce overall groin wound complication rates following peripheral vascular surgery. *J Hosp Infect*. 2018;99:75–80. 3. Nherera LM, Saunders C, Verma S, Trueman P, Fatoye F. Single-use negative pressure wound therapy reduces costs in closed surgical incisions: UK and US economic evaluation. *Journal of Wound Care*. 2021 May 1;30(Sup5):S23–31. 4. Svensson-Björk R, Saha S, Acosta S, Gerdtham UG, Hasselmann J, Asciutto G, Zarrouk M. Cost-effectiveness analysis of negative pressure wound therapy dressings after open inguinal vascular surgery—The randomised INVIPS-Trial. *Journal of Tissue Viability*. 2021 Feb 1;30(1):95–101. 5. Aicher B, Curry P, Croal-Abrahams L, Hao S, Kalsi R, Menon N, Drucker C, Harris D, Toursavadkohl S, Crawford R, Rosenberger S. Infrainguinal wound infections in vascular surgery: an antiquated challenge without a modern solution. *Journal of Vascular Nursing*. 2017 Sep 1;35(3):146–56. 6. Gwilym BL, Dovell G, Dattani N, Ambler GK, Shalhoub J, Forsythe RO, Benson RA, Nandhra S, Preece R, Onida S, Hitchman L. Systematic review and meta-analysis of wound adjuncts for the prevention of groin wound surgical site infection in arterial surgery. *European Journal of Vascular and Endovascular Surgery*. 2021 Jan 7. 7. Grice EA, Segre JA. The skin microbiome. *Nature reviews microbiology*. 2011 Apr;9(4):244–53. 8. Turtiainen J, Hakala T, Hakkarainen T, Karhukorpi J. The impact of surgical wound bacterial colonization on the incidence of surgical site infection after lower limb vascular surgery: a prospective observational study. *European Journal of Vascular and Endovascular Surgery*. 2014 Apr 1;47(4):411–7. 9. Totty JP, Moss JW, Barker E, Mealing SJ, Posnett JW, Chetter IC, Smith GE. The impact of surgical site infection on hospitalisation, treatment costs, and health related quality of life after vascular surgery. *International Wound Journal*. 2021 Jun;18(3):261–8. 10. Audu CO, Columbo JA, Sun SJ, Perri JL, Goodney PP, Stone DH, Nolan BW, Suckow BD. Variation in timing and type of groin wound complications highlights the need for uniform reporting standards. *Journal of vascular surgery*. 2019 Feb 1;69(2):532–43. 11. Siracuse JJ, Nandivada P, Giles KA, Hamdan AD, Wyers MC, Chaikof EL, Pomposelli FB, Schermerhorn ML. Prosthetic graft infections involving the femoral artery. *Journal of vascular surgery*. 2013 Mar 1;57(3):700–5. 12. World Union of Wound Healing Societies (WUWHs) Consensus Document. Closed surgical incision management: understanding the role of NPWT. *Wounds International*. 2016. 13. The World Health Organisation (2016) Global Guidelines for the Prevention of Surgical Site Infection. [online] Available from: <https://apps.who.int/iris/bitstream/handle/10665/250680/9789241549882-eng.pdf?sequence=8> Last Accessed December 2021. 14. Hasselmann J, Björk J, Svensson-Björk R, Acosta S. Inguinal Vascular Surgical Wound Protection by Incisional Negative Pressure Wound Therapy—A Randomized Controlled Trial—INVIPS Trial. *European Journal of Vascular and Endovascular Surgery*. 2019 Dec 1;58(6):e726–7. 15. Wikkeling M, Mans J, Styche T. Single use negative pressure wound therapy in vascular patients: clinical and economic outcomes. *Journal of Wound Care*. 2021 Sep 2;30(9):705–10. 16. Dowsett C, Hampton K, Myers D, Styche T. Use of PICO to improve clinical and economic outcomes in hard-to-heal wounds. *Wounds International*. 2017;8(2):52–58. 17. Kirsner R, Dove C, Reyzelman A, Vayser D, Jaimes H. A Prospective, Randomised, Controlled Clinical Trial on the Efficacy of a single-use Negative Pressure Wound Therapy System, compared to Traditional Negative Pressure Wound Therapy in the Treatment of Chronic Ulcers of the Lower Extremities. *Wound Repair Regen*. 2019;27(5):519–529. 18. Payne C, Edwards D. Application of the Single Use Negative Pressure Wound Therapy Device (PICO) on a Heterogeneous Group of Surgical and Traumatic Wounds. *ePlasty*. 2014;152–166. 19. Rossington A. 2015. A prospective, open, non-comparative, multicentre study to evaluate the functionality and dressing performance of a new negative pressure enhanced dressing (NPED) in acute wounds. CT09/02, May. 20. Smith+Nephew January 2019. Air Leak Tolerance Report: A comparison of PICO v2 (PICO 7 and PICO 14) Devices to PICO v1.6 (PICO) Devices. Internal Report. RD/19/006. 21. Sharpe A, Myers D, Searle R. Using single use negative pressure wound therapy for patients with complicated diabetic foot ulcers: an economic perspective. *Wounds UK*. 2018;14(3):80–137. 22. Test Report. Internal Report. RD/18/133.Smith+Nephew December 2018. Summary of PICO v2 (PICO 7 and PICO 14) Human Factors Testing. Internal Report. RD/18/136. 23. Smith+Nephew January 2019. Acoustic Testing Report: Comparison of PICO v2 (PICO 7 and PICO 14) Devices to PICO v1.6 (PICO) Devices. Internal Report. RD/19/005. 24. Smith+Nephew 2011. The compatibility of ACTICOAT Flex 3 and PICO dressings under NPWT. Internal Report. 1102008. 25. Canonico S, Campitiello F, Della Corte A. Therapeutic possibilities of portable NPWT. *Acta Vulnologica* 10 (2012): 57–64. 26. Hyldig N, Birke-Sorensen H, Kruse M, Vinter C, Joergensen JS, Sorensen JA, Mogensen O, Lamont RF, Bille C. Meta-analysis of negative-pressure wound therapy for closed surgical incisions. *British Journal of Surgery*. 2016 Apr;103(5):477–86. 27. Loveluck J, Copeland T, Hill J, Hunt A, Martin R. Biomechanical modeling of the forces applied to closed incisions during single-use negative pressure wound therapy. *ePlasty*. 2016;16. 28. Malmström M, Huddleston E, Martin R. Biological effects of a disposable, canisterless negative pressure wound therapy system. *ePlasty*. 2014;14. 29. Pellino G, Sciaudone G, Candilio G, Campitiello F, Selvaggi F, Canonico S. Effects of a new pocket device for negative pressure wound therapy on surgical wounds of patients affected with Crohn's disease: a pilot trial. *Surgical Innovation*. 2014 Apr;21(2):204–12. 30. Smith+Nephew 2020. Bacterial barrier testing of the PICO dressing. Internal Report. 2001002. 31. Kilpadi DV, Cunningham MR. Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system. *Wound Repair Regen*. 2011;19(5):588–596. 32. Ma Z, Shou K, Li Z, et al. Negative pressure wound therapy promotes vessel destabilization and maturation at various stages of wound healing and thus influences wound prognosis. *Exp Ther Med*. 2016;11(4):1307–1317. 33. Xia CY, Yu AX, Qi B, et al. Analysis of blood flow and local expression of angiogenesis associated growth factors in infected wounds treated with negative pressure wound therapy. *Mol Med Rep*. 2014;9(5):1749–1754. 34. Birke-Sorensen H, Malmström M, Rome P, et al. Evidence-based recommendations for negative pressure wound therapy: treatment variables (pressure levels, wound filler and contact layer)-steps towards an international consensus. *J Plast Reconstr Aesthet Surg*. 2011;64 Suppl:S1–16. 35. Scalise A, Calamita R, Tartaglione C, et al. Improving wound healing and preventing surgical site complications of closed surgical incisions: a possible role of Incisional Negative Pressure Wound Therapy. A systematic review of the literature. *Int Wound J*. 2016;13(6):1260–1281. 36. Shim HS, Choi JS, Kim SW. A Role for Postoperative Negative Pressure Wound Therapy in Multitissue Hand Injuries. *Biomed Res Int*. 2018. 37. Hudson DA, Adams KG, Van Huyssteen A, Martin R, Huddleston EM. Simplified negative pressure wound therapy: clinical evaluation of an ultraportable, no-cannister system. *Int Wound J*. 2015;12(2):195–201. 38. Stryja J, Staffa R, Říha D, Stryjová K, Nicielniková K. Cost-effectiveness of negative pressure wound therapy in outpatient setting. *Prolekare*. 2015;94(8):322–328. 39. Smith+Nephew 2020. Bacterial barrier testing of the PICO dressing. Internal Report. 2001002. 40. Smith+Nephew 2018. Summary of routine QA testing on MVP of PICO dressings. 2018. Internal Report. DS/18/153/R. 41. McManus H, Woodmansey E. Bacterial retention within a multi-layered absorbent AIRLOCK™ Technology Single Use Negative Pressure Wound Therapy (sNPWT) dressing. Paper presented at: EWMA; 2018; Krakow, Poland. 42. Smith+Nephew July 2018. PICO 7Y Non-NPWT Wound Model Summary. Internal Report. DS.18.260.R. 43. Casey C. Consistent delivery of therapeutic negative pressure levels by a single use negative pressure wound therapy system (sNPWT)* in a wound model. Paper presented at: EWMA; 2019; Gothenburg, Sweden. 44. Smith+Nephew 2021. PICO® Pressure Mapping Study. Internal Report. DS/19/211/R - Part B. 45. NICE (2019) PICO negative pressure wound dressings for closed surgical incisions [online] accessible from: <https://www.nice.org.uk/guidance/mtg43> Last accessed December 2021. 46. Smith+Nephew 2019. PICO Biomechanical Study. Internal Report. DS/19/211/R.



HOW ANATOMY DICTATES OUTCOMES



BACK AT HOME, NOT IN HOSPITAL



THE PICO SYSTEM



ORDERING




REFERENCES


PICO® sNPWT FOR VASCULAR SURGERY | 10

Negative pressure wound therapy (NPWT)

NPWT has **multiple mechanisms of action** which may help promote incisional wound healing and **reduce the odds of SSCs**.²⁵⁻³⁰



*As demonstrated in biomechanical modelling †As demonstrated in vivo



The PICO[◇] System uses AIRLOCK[◇] Technology

AIRLOCK Technology **ensures consistent delivery** of negative pressure,* protecting the incision and treating the wider zone of injury.^{17,30,41-44,46}
Only PICO sNPWT dressings have AIRLOCK Technology.



Wide delivery.
Constant pressure.
Optimal outcomes.

Super absorbent core locking exudate away from wound⁴²

Top film layer has a high moisture vapour transmission rate and protects the wounds from external contamination^{39,40}

Soft Port with integrated filter

Silicone adhesive layer helps to minimise pain on removal^{18,37,38}

Pioneering **AIRLOCK Technology** transmits pressure across a wider zone of treatment beyond the incision itself^{43,44,46}

Up to
80%
of the exudate is lost by evaporation²⁸

While
20%
is absorbed in the dressing²⁸

*as demonstrated in benchtop testing

High quality evidence for high risk patients

In a meta-analysis¹ of **29 studies**, including **11 randomised controlled trials** (RCTs) with a total of **5,614 patients**, in a variety of surgical incisions following orthopaedics, vascular, breast and obstetrics surgery, the results showed:



For the NHS in the UK, the PICO System was found to provide clinical benefit and **save £105 per patient**.^{||3}
Conversion \$162.47CDN*



Recent NICE guidance demonstrates that PICO sNPWT provides better outcomes than standard care for preventing surgical site complications in high-risk patients with closed surgical incisions.⁴⁵

Incremental acquisition costs of PICO sNPWT are more than offset by savings in the treatment of SSIs.⁴⁵

*Compared to care with standard dressings; p<0.025; meta-analysis of 19 studies (odds ratio (OR): 0.37)
†Compared to care with standard dressings; p<0.025; meta-analysis of 9 studies (OR: 0.70)
‡Compared to care with standard dressings; p<0.025; meta-analysis of 6 studies (OR: 0.23)
§Compared to care with standard dressings; p<0.025; meta-analysis of 10 studies
||over a 12-week period compared with standard care
*Conversion CDN July 2022