

Introducing the new **RENASYS[◇] EDGE** Negative Pressure Wound Therapy device

Easy to use:

- Pump offers an intuitive, user-friendly interface for easy to learn operation and troubleshooting¹

Built for efficiency:

- Sturdy, reliable and durable to help minimise the number of device returns²

Patient-friendly:

- Size and weight of device is designed to allow patients to continue with their daily lives and activities³

Smith+Nephew

RENASYS[◇] EDGE

Negative Pressure Wound
Therapy System

smith-nephew.com/renasysedge



The RENASYS[◇] EDGE System



RENASYS[◇] EDGE
Negative Pressure Wound
Therapy System

Comparing the systems:

GO vs TOUCH vs EDGE

The **RENASYS[®] EDGE** pump:

- A smaller footprint
- A longer battery life
- Less maintenance

RENASYS GO System



RENASYS TOUCH System






RENASYS EDGE System





	GO	TOUCH	EDGE
Weight	1.13 Kg	1.1Kg	1.13 Kg
Dimensions	17.5cmx24cm	18cmx19cm	15cmx11cm
Footprint area	420 cm ²	N/A	165 cm ²
Need annual service and maintenance	Yes	Yes	No
Canister with solidifier	300ml or 750ml	300ml or 800ml	300ml hidden or 800ml
Integrated carry handle	No	No	Yes
Self-standing	No	Kickstand with 300ml canister	Yes
Display	Monochromatic screen	Full-colour screen	Full-color screen
Operating time	Up to 20 hours (therapy)	10-16 hours (therapy)	> 24 hours (therapy)
Pressure range	-40 to -200mmHg	-25 to -200mmHg	-25 to -200mmHg
Therapy settings	Continuous and Intermittent	Continuous and intermittent	Continuous and Intermittent
Alarm volumes	Fixed 46dB	Selectable High: 74dB Medium: 68dB Low: 60dB	Selectable High: 75dB Medium: 65dB Low: 55dB
Adjustable draw down rates while therapy starts	No	Yes	Yes
On-board tutorials	No	Yes	Yes
Pump self-test functionality	No	No	Yes
O-ring	Serviceable part. Integrated in the pump	Serviceable part. Integrated in the pump.	Non-serviceable part. Integrated in disposable canister
Odour filter	Serviceable part. Integrated in the pump	Serviceable part. Integrated in the pump.	Non-serviceable part. Integrated in disposable canister
Home healthcare lock	Patient mode restricts access to therapy settings	Patient mode restricts access to therapy settings	Pump locks automatically
Autostart after therapy paused for too long	No	No	Yes – after 1 hour of inactivity
Accessories	Carry bag and carry straps	Carry bag and carry straps	Carry bag and carry straps
Pump NFC sensor linked to on-line resources	No	No	Yes

*Compared to RENASYS TOUCH and GO Systems

Ordering codes:

Item description	Item SKU	
Pump	66803126	
300ml canister with solidifier	66803139	
800ml canister with solidifier	66803140	

Item description	Item SKU	
Carrying bag	66803136	
Transit case	66803138	
Clinician manual	66803143	
Homecare manual	66803144	

Advanced Wound Management
 Smith & Nephew, Inc.
 2280 Argentia Rd.
 Mississauga, ON
 L5N 6H8, CAN
 T 1-800-463-7439
 F 1-800-671-9140
 NPWT 24/7 Clinical Hotline:
 1-800-463-7439

www.smith-nephew.com

♦Trademark of Smith+Nephew
 All Trademarks acknowledged
 © March 2023 Smith+Nephew
 CA99099 04/23 | GMC1684

Shaping what's possible in wound care

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

References: **1.** Smith+Nephew 2022. RENASYS EDGE System Human Factors Summative Report Summary. Internal Report. CSD.AWM.22.071. **2.** Smith+Nephew 2022. Internal Report. CSD.AWM.22.069. **3.** Smith+Nephew 2022. Summary of footprint, portability, wearability, weight and audible noise for the RENASYS[®] EDGE system. Internal Report. CSD.AWM.22.067. **4.** Smith+Nephew 2022. How the RENASYS EDGE Negative Pressure Wound Therapy System provides continuity of care to the patient Internal Report. EO.AWM.PCS270.003.V1.