


Advance Serenity[®] 14

HYDROPHILIC PTA BALLOON DILATATION CATHETER

A smoother crossing¹



Next-gen hydrophilic coating

Provides lubricity to navigate with low friction.



Laser-formed tip

Creates trackability and smooth transitions.¹



Coaxial shaft design

Balances stiffness and flexibility to deliver pushability.¹



Pleated balloon

Facilitates a smooth rewrap.



Order Number	Reference Part Number	Wire Guide Diameter inch	Balloon Inflated Diameter mm	Balloon Length cm	Recommended Introducer Fr	Exact Catheter Diameter Fr	Nominal Pressure atm	Rated Burst Pressure atm
0.014 inch Diameter Guide – 120 cm Long OTW Shaft								
G57848	SER14-120-20-20	0.014	2	2	4	3.9	6	16
G57850	SER14-120-20-40	0.014	2	4	4	3.9	6	16
G57849	SER14-120-20-60	0.014	2	6	4	3.9	6	16
G57851	SER14-120-20-80	0.014	2	8	4	3.9	6	16
G57852	SER14-120-20-100	0.014	2	10	4	3.9	6	16
G57853	SER14-120-20-150	0.014	2	15	4	3.9	6	16
G57854	SER14-120-20-220	0.014	2	22	4	3.9	6	16
G57855	SER14-120-25-20	0.014	2.5	2	4	3.9	6	16
G57856	SER14-120-25-40	0.014	2.5	4	4	3.9	6	16
G57857	SER14-120-25-60	0.014	2.5	6	4	3.9	6	16
G57858	SER14-120-25-80	0.014	2.5	8	4	3.9	6	16
G57859	SER14-120-25-100	0.014	2.5	10	4	3.9	6	16
G57860	SER14-120-25-150	0.014	2.5	15	4	3.9	6	16
G57861	SER14-120-25-220	0.014	2.5	22	4	3.9	6	16
G57862	SER14-120-30-20	0.014	3	2	4	3.9	6	16
G57863	SER14-120-30-40	0.014	3	4	4	3.9	6	16
G57864	SER14-120-30-60	0.014	3	6	4	3.9	6	16
G57865	SER14-120-30-80	0.014	3	8	4	3.9	6	16
G57866	SER14-120-30-100	0.014	3	10	4	3.9	6	16
G57867	SER14-120-30-150	0.014	3	15	4	3.9	6	16
G57868	SER14-120-30-220	0.014	3	22	4	3.9	6	16
G57869	SER14-120-35-20	0.014	3.5	2	4	3.9	6	16
G57870	SER14-120-35-40	0.014	3.5	4	4	3.9	6	16
G57871	SER14-120-35-60	0.014	3.5	6	4	3.9	6	16
G57872	SER14-120-35-80	0.014	3.5	8	4	3.9	6	16
G57873	SER14-120-35-100	0.014	3.5	10	4	3.9	6	16
G57874	SER14-120-35-150	0.014	3.5	15	4	3.9	6	16
G57875	SER14-120-35-220	0.014	3.5	22	4	3.9	6	16
G57876	SER14-120-40-20	0.014	4	2	4	3.9	6	16
G57877	SER14-120-40-40	0.014	4	4	4	3.9	6	16
G57878	SER14-120-40-60	0.014	4	6	4	3.9	6	16
G57879	SER14-120-40-80	0.014	4	8	4	3.9	6	16
G57880	SER14-120-40-100	0.014	4	10	4	3.9	6	16
G57881	SER14-120-40-150	0.014	4	15	4	3.9	6	16
G57882	SER14-120-40-220	0.014	4	22	4	3.9	6	16

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

Order Number	Reference Part Number	Wire Guide Diameter inch	Balloon Inflated Diameter mm	Balloon Length cm	Recommended Introducer Fr	Exact Catheter Diameter Fr	Nominal Pressure atm	Rated Burst Pressure atm
0.014 inch Diameter Guide – 150 cm Long OTW Shaft								
G57883	SER14-150-20-20	0.014	2	2	4	3.9	6	16
G57884	SER14-150-20-40	0.014	2	4	4	3.9	6	16
G57885	SER14-150-20-60	0.014	2	6	4	3.9	6	16
G57886	SER14-150-20-80	0.014	2	8	4	3.9	6	16
G57887	SER14-150-20-100	0.014	2	10	4	3.9	6	16
G57888	SER14-150-20-150	0.014	2	15	4	3.9	6	16
G57889	SER14-150-20-220	0.014	2	22	4	3.9	6	16
G57890	SER14-150-25-20	0.014	2.5	2	4	3.9	6	16
G57891	SER14-150-25-40	0.014	2.5	4	4	3.9	6	16
G57892	SER14-150-25-60	0.014	2.5	6	4	3.9	6	16
G57893	SER14-150-25-80	0.014	2.5	8	4	3.9	6	16
G57894	SER14-150-25-100	0.014	2.5	10	4	3.9	6	16
G57895	SER14-150-25-150	0.014	2.5	15	4	3.9	6	16
G57896	SER14-150-25-220	0.014	2.5	22	4	3.9	6	16
G57897	SER14-150-30-20	0.014	3	2	4	3.9	6	16
G57898	SER14-150-30-40	0.014	3	4	4	3.9	6	16
G57899	SER14-150-30-60	0.014	3	6	4	3.9	6	16
G57900	SER14-150-30-80	0.014	3	8	4	3.9	6	16
G57901	SER14-150-30-100	0.014	3	10	4	3.9	6	16
G57902	SER14-150-30-150	0.014	3	15	4	3.9	6	16
G57903	SER14-150-30-220	0.014	3	22	4	3.9	6	16
G57904	SER14-150-35-20	0.014	3.5	2	4	3.9	6	16
G57905	SER14-150-35-40	0.014	3.5	4	4	3.9	6	16
G57906	SER14-150-35-60	0.014	3.5	6	4	3.9	6	16
G57907	SER14-150-35-80	0.014	3.5	8	4	3.9	6	16
G57908	SER14-150-35-100	0.014	3.5	10	4	3.9	6	16
G57909	SER14-150-35-150	0.014	3.5	15	4	3.9	6	16
G57910	SER14-150-35-220	0.014	3.5	22	4	3.9	6	16
G57911	SER14-150-40-20	0.014	4	2	4	3.9	6	16
G57912	SER14-150-40-40	0.014	4	4	4	3.9	6	16
G57913	SER14-150-40-60	0.014	4	6	4	3.9	6	16
G57914	SER14-150-40-80	0.014	4	8	4	3.9	6	16
G57915	SER14-150-40-100	0.014	4	10	4	3.9	6	16
G57916	SER14-150-40-150	0.014	4	15	4	3.9	6	16
G57917	SER14-150-40-220	0.014	4	22	4	3.9	6	16

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.



Advance Serenity® 14


HYDROPHILIC PTA BALLOON DILATATION CATHETER

The Advance Serenity 14 Hydrophilic PTA Balloon Dilatation Catheter is indicated for percutaneous transluminal angioplasty (PTA) dilation of peripheral vasculature stenosis in the iliac, femoral, ilio-femoral, popliteal, infrapopliteal, and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae.

- All of the balloon sizes are compatible with 4 Fr sheaths.
- Balloon lengths of up to 220 mm allow for the treatment of longer lesions with fewer inflations and deflations.

Benchtop test crossing 0.4 inch synthetic lesion

Advance Serenity 14 vs other 0.014 inch PTA balloon catheters²



2. Advance Serenity. Test Report RD2016-09. (25 APR 2016). Surmodics, Inc.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 8007615, IS.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-8.5X11