


AZUR[®]

Peripheral HydroCoil Embolization System

AZUR[®] **CX**

Peripheral Coil System



INCREASE YOUR
**EMBOLIZATION
EFFICIENCY**

Use fewer coils and achieve cost savings
with patented Hydrogel technology¹



**TERUMO**
INTERVENTIONAL
SYSTEMS

PUSHING BOUNDARIES

Terumo Interventional Systems is **committed to your success** with innovative procedural solutions and ongoing support for your most challenging cases.

We are relentlessly seeking new ways to help you apply effective solutions and achieve **better outcomes for more patients.**



TERUMO
INTERVENTIONAL
SYSTEMS

AZUR[®]

Peripheral HydroCoil Embolization System

AZUR[®] CX

Peripheral Coil System

ACHIEVE GREATER EFFICIENCIES

Reach your embolization endpoints with superior filling volume and packing density²⁻⁴

Coil volume per 1 cm of length*

AZUR[®] HydroCoil 0.035"



AZUR[®] HydroCoil 0.018"



AZUR[®] CX 0.035"



Boston Scientific Interlock[®] Coil 0.035"



Penumbra Ruby[®] Coil



Medtronic Concerto PGLA[®] Coil




*Data from AngioCalc.com

EMBOLIZE WITH VOLUME IN MIND

Coil length alone does not determine filling volume and packing density

The AZUR® HydroCoil System provides almost 3 times the volume in half the coil length of a bare metal coil.*



*Data from AngioCalc.com

ACHIEVE FILLING VOLUME WITH **FEWER COILS**

Realize cost savings vs. competitive coils

Number of coils and total length needed for a 2 cm spherical aneurysm*

AZUR® HydroCoil 0.035"



2.9 coils
87 cm length

AZUR®
Peripheral HydroCoil Embolization System

AZUR® HydroCoil 0.018"



5.7 coils
172 cm length

AZUR®
Peripheral HydroCoil Embolization System

AZUR® CX 0.035"



6.1 coils
236 cm length

AZUR® CX
Peripheral Coil System

Penumbra Ruby® Coil



8.3 coils
496 cm length

Boston Scientific Interlock® Coil 0.035"



11.3 coils
450 cm length

Medtronic Concerto PGLA® Coil



18.9 coils
944 cm length

*Data from AngioCalc.com

Total volume: 4,189 mm³; Filling requirement 24% (1,055 mm³)

AZUR[®]
Peripheral HydroCoil Embolization System

AZUR[®] CX
Peripheral Coil System

Increase your procedural efficiencies with
a comprehensive line of peripheral coils



AZUR[®] Framing Coil

Available in 0.035" and 0.018"
Detachable Systems



AZUR[®] HydroCoil

Available in 0.035" and 0.018"
Detachable and Pushable Systems



AZUR[®] CX

Available in 0.035" and 0.018"
Detachable Systems

FIND OUT MORE



US: 800.862.4143
CA: 833.883.7866



US: terumo.com
CA: terumocanada.ca

References:

1. Data on File. Terumo Corporation. IMS Data. **2.** Ding YH, Dai D, Lewis DA, Cloft HJ, Kallmes DF. Angiographic and histologic analysis of experimental aneurysms embolized with platinum coils, Matrix, and HydroCoil. *AJNR Am J Neuroradiol.* 2005 Aug;26(7):1757-63. (animal study). **3.** López-Benítez, R., Hallscheidt, P., Kratochwil, C. et al. *Cardiovasc Intervent Radiol* (2013) 36: 105. **4.** Maleux G, Deroose C, Fieuwis S, Van Cutsem E, Heye S, Bosmans H, Verslype C. Prospective comparison of hydrogel-coated microcoils versus fibered platinum microcoils in the prophylactic embolization of the gastroduodenal artery before yttrium-90 radioembolization. *J Vasc Interv Radiol.* 2013 Jun;24(6):797-803.

RX ONLY. The advertisement is directed to physicians only, and not to consumers. Refer to product labels and packaging insert for complete warnings, precautions, potential complications, and instructions for use. Products may not have regulatory approval in all countries. Please contact your local sales representative if you have questions about the availability of products in your area.

©2019 Terumo Medical Corporation and Terumo Medical Canada Inc.
All rights reserved. All brand names are trademarks of Terumo. PM-02478

TERUMO
INTERVENTIONAL
SYSTEMS